

DRI<sup>®</sup> Cannabinoid Assay**Kit Configuration**

P/N W150185	1 x 100 mL Antibody/Substrate Reagent A (R1)
	1 x 100 mL Enzyme Conjugate Reagent E (R2)

**Reagent Preparation**

P/N W150185: Reagents are ready to use. Pour R1 and R2 in the appropriate bottles and place them in the reagent tray.

**In use Stability**

On Board: 30 days

**Specimen**

Urine

**Calibration**

Use: Negative Calibrator Cat. No.W151664

THC 20 ng/mL Calibrator Cat. No.W150235

THC 50 ng/mL Calibrator Cat. No.W150042

THC 100 ng/mL Calibrator Cat. No.W150044

THC 200 ng/mL Calibrator Cat. No.W150206

Recalibrate every 30 days or when a new lot of reagent is used.

**Quality Control**

THC 40 ng/mL Control Cat. No.W150170

THC 60 ng/mL Control Cat. No.W150168

THC 75 ng/mL Control Cat. No.W150214

THC 125 ng/mL Control Cat. No.W150212

**Calculation and Analytical Results**

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the ILab 350 Operator Manual.

**Semiquantitative results**

A rough estimate of drug concentration in the samples can be obtained by running a standard curve with all calibrators and quantifying samples off the standard curve.

**References / Literatur / Bibliografía / Bibliographie / Bibliografia /**

See package insert inclosed in the kit

<b>Performance Characteristics</b>
------------------------------------

The performance below were obtained working with a cutoff of 50 ng/mL

**Limitation/Interfering Substances**

A positive result by this assay should be confirmed by another nonimmunological method such as GC, TLC or GC/MS.

It is possible that other substances and/or factors (eg, technical or procedural) not listed in the specificity table (see package insert) may interfere with the test and cause false results.

**Precision**

	Samples/Runs	Mean (ng/mL)	CV(%)	Mean (ng/mL)	CV(%)	Mean (ng/mL)	CV%
Within Run	5/10	32	5.6	45	4.7	63	2.0
Total	5/10	32	9.5	49	6.9	63	2.4

**Minimum Detection Limit**

1 ng/mL



### Instrument Settings

<b>Chemistry Parameters</b>		R1	
Method	<input type="text"/>	Reagent Name	<input type="text" value="THC"/>
Name	<input type="text" value="THC"/>	Volume	<input type="text" value="100 μL"/>
Unit	<input type="text" value="ng/mL"/>	R2	<input type="text" value="enable"/>
Assay Type	<input type="text" value="Rate"/>	Reagent Name	<input type="text" value="THC"/>
		Volume	<input type="text" value="100 μL"/>
		Wash	<input type="text" value="disable"/>
		Reagent Name	<input type="text"/>
		Reagent Type	<input type="text"/>
		Diluent	<input type="text" value="disable"/>
		Reagent Name	<input type="text"/>
Measuring Points	1 <input type="text" value="disable"/>	start	<input type="text"/>
		end	<input type="text"/>
	2 <input type="text" value="enable"/>	start	<input type="text" value="15"/>
		end	<input type="text" value="18"/>
Wave Length		Decimal Points	<input type="text" value="0"/>
Prim	<input type="text" value="340"/>	Sec	<input type="text" value="405"/>
		Technical Range (Conc)	<input type="text" value="0 - 200"/>
		mAbs/10	<input type="text" value="-30000 / 30000"/>
Sampling Volume	<input type="text" value="10 μL"/>	RPT Wash (R1)	<input type="text" value="Sys Water"/>
Dilution	<input type="text" value="disable"/>	(R2)	<input type="text" value="Sys Water"/>
Rerun ( High)	<input type="text"/>	Instrument Factor a	<input type="text" value="1"/>
Dilution	<input type="text" value="disable"/>	Stirring Speed	R1 <input type="text" value="mid"/> R2 <input type="text" value="0"/>
Rerun ( Low)	<input type="text"/>		

### Calibration Checks

** Duplicate Limit	<input type="text"/>	** mAbs/10	Sampling Method for Standards
** Activity Limit	<input type="text"/>	** mAbs/10	<input type="text"/> Duplicate
			<input checked="" type="checkbox"/> Triplicate
** Linearity Limit	<input type="text"/>	** %	Blank measurement
** Prozone Limit	<input type="text" value="upper"/>		<input checked="" type="checkbox"/> Disable Reagent blank
SL1-S	** <input type="text" value="SL1-F"/>	**	<input type="text" value="None"/>
SL2-S	** <input type="text" value="SL2-F"/>	**	Reagent blank measurement at calibration
Sens	<input type="text"/>	mAbs/10	<input type="text"/> Reagent blank (system water)
<input checked="" type="checkbox"/> Absorbance Limit			** Multiplex measurement is the same as standards
Reaction	<input type="text" value="Increase"/>		Reagent Blank Limit Checks
Limit	<input type="text" value="25000"/>	mAbs/10	** Duplicate limit <input type="text" value="50"/> mAbs/10

### Calibration

Method	<input type="text"/>	Name	<input type="text" value="THC"/>	Interval	<input type="text" value="30"/> days
Calculation	<input type="text" value="Point to Point"/>				
	Conc	WORK	MASTER	Lot No	
S1	<input type="text" value="0"/>				K <input type="text" value="N/A"/>
S2	<input type="text" value="20"/>				
S3	<input type="text" value="50"/>				
S4	<input type="text" value="100"/>				
S5	<input type="text" value="200"/>				
S6	<input type="text"/>				

### Reagent Registration

Reagent Code	<input type="text" value="0146"/>				
Reagent Name	<input type="text" value="THC"/>				
R1	<input checked="" type="checkbox"/> enable	Volume (L)	<input type="text" value="**"/> mL	Volume (S)	<input type="text" value="**"/> mL
R2	<input checked="" type="checkbox"/> enable	Volume (L)	<input type="text" value="**"/> mL	Volume (S)	<input type="text" value="**"/> mL
		Stability Check	<input checked="" type="checkbox"/> enable	Term	<input type="text" value="30"/> days
		Stability Check	<input checked="" type="checkbox"/> enable	Term	<input type="text" value="30"/> days

\*\* Operator definable N/A not applicable to this test