

Kit Configuration

P/N 3000-2238	2 x 16 mL IgE R1
	2 x 5 mL IgE R2

Reagent Preparation

P/N 3000-2238	IgE R1: Ready to use
	IgE R2: Ready to use. Invert to mix well before first use. Avoid foam formation.
	Place the bottles into reagent tray.

In Use Stability

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial securely closed.

Specimen

Serum.

Calibration

Use quantex IgE standard multipoint Cat. No 300-2240. See vial labels for lot specific concentrations. A reagent blank should be run daily before sample analysis. Recalibrate every 72 days or when a new lot of reagent is used.

Quality Control

Use quantex Ferritin/Myoglobin/IgE Control I/II Cat. No. 3000-2222.

Calculation of Analytical Results

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the ILab 350 Operator Manual.

Reference Interval

Age related concentrations must be taken into account when interpreting IgE values in children. IgE does not cross the placental barrier so IgE is not detectable in new borns. The IgE concentration increases during first years of life, reaching a peak at 10-15 years and dropping subsequently to adult values.

Age	1-12 months	1-5 years	6-9years	10-15years	Adults
Concentration (IU/mL)	<15	<60	<90	<200	<100

In any case, these concentrations are only indicatives and each laboratory should establish its own reference range.

References / Literatur / Bibliografia / Bibliographie / Bibliografia /

See package insert enclosed in the kit

Performance Characteristics

Limitation/Interfering Substances

No significant interference from bilirubin up to concentrations of 14.7 mg/dL (250 µmol/L) and hemoglobin up to concentrations of 1600 mg/dL (0.96 mmol/L). For a comprehensive review of interfering substances, refer to the publication by Young *et al.*¹

Precision

	Samples/Runs	Mean (IU/mL)	CV(%)	Mean (IU/mL)	CV(%)
Within run	4/10	66	7.3	485	1.0
Total	4/10	66	11.9	485	3.0

Method Comparison

Comparison Method (x)	same reagent
Comparison Instrument (x)	Nephelometer (Behring BN II)
Slope	0.893
y intercept	-0.264
Mean X (IU/mL)	175
Mean Y (IU/mL)	156
r	0.994
n	102

Linearity

no rerun 25 -1000 IU/mL ; with rerun 25 - 1000 IU/mL

Minimum Detection Limit

18 IU/mL

Quantification Limit

25 IU/mL

Instrument Settings

Chemistry Parameters				R1	
Method	<input type="text"/>	Reagent Name	<input type="text" value="IgE"/>	Volume	<input type="text" value="200 μL"/>
Name	<input type="text" value="IgE"/>	R2	<input type="text" value="enable"/>		
Unit	<input type="text" value="IU/mL"/>	Reagent Name	<input type="text" value="IgE"/>	Volume	<input type="text" value="80 μL"/>
Assay Type	<input type="text" value="End"/>	Wash	<input type="text" value="disable"/>	Reagent Name	<input type="text"/>
		Diluent	<input type="text" value="enable"/>	Reagent Type	<input type="text"/>
				Reagent Name	<input type="text" value="Saline"/>
Measuring Points	1 enable	start	<input type="text" value="14"/>	Decimal Points	<input type="text" value="0"/>
		end	<input type="text" value="14"/>		
	2 enable	start	<input type="text" value="21"/>	Normal Range	<input type="text" value="0"/> <input type="text" value="200"/>
		end	<input type="text" value="21"/>		
Wave Length				Technical Range (Conc)	<input type="text" value="0.0"/> <input type="text" value="1000"/>
Prim	<input type="text" value="570"/>	Sec	<input type="text"/>	mAbs/10	<input type="text" value="-30000"/> <input type="text" value="30000"/>
Sampling Volume	<input type="text" value="5 μL"/>				
Dilution	<input type="text" value="disable"/>	RPT Wash	(R1) <input type="text" value="Sys Water"/>		
	<input type="text" value="μL"/> <input type="text" value="μL"/>		(R2) <input type="text" value="Sys Water"/>		
Rerun (High)	<input type="text" value="5 μL"/>				
Dilution	<input type="text" value="enable"/>	Instrument Factor a	<input type="text" value="1"/>	b	<input type="text" value="0"/>
	<input type="text" value="15 μL"/> <input type="text" value="135 μL"/>	Stirring Speed	R1 <input type="text" value="high"/>	R2	<input type="text" value="high"/>
Rerun (Low)	<input type="text" value="25 μL"/>				

Calibration Checks

** Duplicate Limit	<input type="text"/>	** mAbs/10	Sampling Method for Standards	
** Sensitivity Limit	<input type="text"/>	** mAbs/10	<input checked="" type="checkbox"/> Duplicate	
			<input type="checkbox"/> Triplicate	
** Linearity Limit	<input type="text"/>	** %	Blank measurement	
** Prozone Limit	<input type="text" value="upper"/>		<input checked="" type="checkbox"/> Enable Reagent blank	
SL1-S	** <input type="text" value="SL1-F"/>	**	<input type="text" value="None"/>	
SL2-S	** <input type="text" value="SL2-F"/>	**	Reagent blank measurement at calibration	
Sens	<input type="text"/>	mAbs/10	<input checked="" type="checkbox"/> Reagent blank (system water)	
<input checked="" type="checkbox"/> Absorbance Limit			Multiplex measurement is the same as standards	
Reaction	<input type="text" value="Increase"/>		Reagent Blank Limit Checks	
Limit	<input type="text" value="25000"/>	mAbs/10	** Duplicate limit	<input type="text" value="50"/> mAbs/10

Calibration

Method	<input type="text"/>	Name	<input type="text" value="IgE"/>	Interval	<input type="text" value="72"/> days
Calculation	<input type="text" value="Point to Point"/>				
	Conc	WORK	MASTER	Lot No	
S1	<input type="text" value="0"/>	<input type="text" value="-10"/>			K <input type="text" value="N/A"/>
S2	<input type="text" value="50"/>	<input type="text" value="132"/>			
S3	<input type="text" value="100"/>	<input type="text" value="270"/>			
S4	<input type="text" value="200"/>	<input type="text" value="599"/>			
S5	<input type="text" value="500"/>	<input type="text" value="1983"/>			
S6	<input type="text" value="1000"/>	<input type="text" value="4639"/>			

Reagent Registration

Reagent Code	<input type="text" value="0189"/>				
Reagent Name	<input type="text" value="IgE"/>				
R1	<input checked="" type="checkbox"/> enable	Volume (L)	<input type="text" value="**"/> mL	Volume (S)	<input type="text" value="**"/> mL
R2	<input checked="" type="checkbox"/> enable	Volume (L)	<input type="text" value="**"/> mL	Volume (S)	<input type="text" value="**"/> mL
		Stability Check	<input checked="" type="checkbox"/> enable	Term	<input type="text" value="**"/> days
		Stability Check	<input checked="" type="checkbox"/> enable	Term	<input type="text" value="**"/> days
**	Operator definable	N/A	not applicable to this test	Calibration curve is only as example	