

Kit Configuration

| | |
|---------------|-----------------|
| P/N 3000-2255 | 1 x 75 mL RF R1 |
| | 2 x 6 mL RF R2 |

Reagent Preparation

| | |
|---------------|---------------------------------------------------------------------------------|
| P/N 3000-2255 | RF R1: Ready to use. |
| | RF R2: Ready to use. Invert to mix well before first use. Avoid foam formation. |
| | Place the bottles into reagent tray. |

In Use Stability

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial securely closed.

Specimen

Serum.

Calibration

Use quantex RF Plus standard Cat. No 300-2252. The calibrator contains 200 IU/mL of RF. Serial dilution of standard with saline automatically by instrument as follows: 0 IU/mL (saline) 25, 50, 100 and 200 IU/mL (undiluted standard). A reagent blank should be run daily before sample analysis. Recalibrate every 90 days or when a new lot of reagent is used.

Quality Control

Use quantex ASO/CRP/RF Control I Cat. No. 3000-2069 and quantex ASO/CRP/RF Control II Cat. No. 3000-2070. .

Calculation of Analytical Results

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the ILab 350 Operator Manual.

Reference Interval

Results of RF lower than 30 IU/mL are considered normal. Results between 30 to 50 IU/mL are considered weakly positive.

References / Literatur / Bibliografía / Bibliographie / Bibliografia /

See package insert enclosed in the kit

Performance Characteristics

Limitation/Interfering Substances

No significant interference from lipemia up to a sample absorbance of 1.0/cm at 660 nm or 300 mg/dL triglycerides (3.39 mmol/L) bilirubin up to concentrations of 33 mg/dL (560 µmol/L) , hemoglobin up to concentrations of 150 mg/dL (0.09 mmol/L). For a comprehensive review of interfering substances, refer to the publication by Young *et al.*¹

Precision

| | Samples/Runs | Mean (IU/mL) | CV(%) | Mean (IU/mL) | CV(%) |
|------------|--------------|--------------|-------|--------------|-------|
| Within run | 4/10 | 50 | 1.2 | 131 | 2.3 |
| Total | 4/10 | 50 | 2.6 | 131 | 3.1 |

Method Comparison

| | |
|---------------------------|-----------------|
| Comparison Method (x) | Quantex RF Plus |
| Comparison Instrument (x) | Cobas Mira |
| Slope | 1.00 |
| y intercept | 0.0629 |
| Mean X (IU/mL) | 172 |
| Mean Y (IU/mL) | 171 |
| r | 0.99 |
| n | 40 |

Linearity

no rerun 21 - 200 IU/mL ; with rerun 21 - 2000 IU/mL

Minimum Detection Limit

2 IU/mL

Quantification Limit

21 IU/mL

Instrument Settings

| | | | |
|-----------------------------|------------------------------------------------------------------------------------|------------------------|----------------------------------------------------------------------------------------------------------------------------|
| Chemistry Parameters | | R1 | |
| Method Name | <input type="text" value="RF"/> | Reagent Name | <input type="text" value="RF"/> Volume <input type="text" value="190 μL"/> |
| Unit | <input type="text" value="IU/mL"/> | R2 | <input type="text" value="enable"/> Reagent Name <input type="text" value="RF"/> Volume <input type="text" value="40 μL"/> |
| Assay Type | <input type="text" value="End"/> | Wash | <input type="text" value="disable"/> Reagent Name <input type="text"/> |
| | | Diluent | <input type="text" value="enable"/> Reagent Type <input type="text"/> |
| | | | Reagent Name <input type="text" value="Saline"/> |
| Measuring Points | 1 enable start <input type="text" value="12"/> end <input type="text" value="13"/> | Decimal Points | <input type="text" value="1"/> |
| | 2 enable start <input type="text" value="25"/> end <input type="text" value="26"/> | Normal Range | <input type="text" value="0"/> <input type="text" value="30"/> |
| Wave Length | Prim <input type="text" value="570"/> Sec <input type="text"/> | Technical Range (Conc) | <input type="text" value="0.0"/> <input type="text" value="200"/> |
| | | mAbs/10 | <input type="text" value="-30000 / 30000"/> |
| Sampling Volume | <input type="text" value="6 μL"/> | RPT Wash (R1) | <input type="text" value="Sys Water"/> |
| Dilution | <input type="text" value="disable"/> | (R2) | <input type="text" value="Sys Water"/> |
| Rerun (High) | <input type="text" value="6 μL"/> | Instrument Factor a | <input type="text" value="1"/> <input type="text" value="b 0"/> |
| Dilution | <input type="text" value="enable"/> | Stirring Speed | R1 <input type="text" value="high"/> R2 <input type="text" value="high"/> |
| Rerun (Low) | <input type="text" value="15 μL"/> <input type="text" value="135 μL"/> | | |
| | <input type="text" value="12 μL"/> | | |

Calibration Checks

| | | | |
|------------------------------------------------------|----------------------------------------------------------------------------------------------------|--------------------------------------------------|------------------------------------------------------------------|
| ** Duplicate Limit | <input type="text" value="**"/> | mAbs/10 | Sampling Method for Standards |
| ** Sensitivity Limit | <input type="text" value="**"/> | mAbs/10 | <input checked="" type="checkbox"/> Duplicate |
| | | | <input type="checkbox"/> Triplicate |
| ** Linearity Limit | <input type="text" value="**"/> | % | Blank measurement |
| ** Prozone Limit | <input type="text" value="upper"/> | | <input checked="" type="checkbox"/> Enable Reagent blank |
| SL1-S | <input type="text" value="**"/> <input type="text" value="SL1-F"/> <input type="text" value="**"/> | | <input type="text" value="None"/> |
| SL2-S | <input type="text" value="**"/> <input type="text" value="SL2-F"/> <input type="text" value="**"/> | | Reagent blank measurement at calibration |
| Sens | <input type="text" value="**"/> | mAbs/10 | <input checked="" type="checkbox"/> Reagent blank (system water) |
| <input checked="" type="checkbox"/> Absorbance Limit | Reaction <input type="text" value="Increase"/> | Limit <input type="text" value="25000"/> mAbs/10 | ** Multiplex measurement is the same as standards |
| | | | Reagent Blank Limit Checks |
| | | | ** Duplicate limit <input type="text" value="50"/> mAbs/10 |

Calibration

| | | | | | |
|-------------|---------------------------------------------|------------------------------------|---------------------------------|----------|--------------------------------------|
| Method | <input type="text"/> | Name | <input type="text" value="RF"/> | Interval | <input type="text" value="90"/> days |
| Calculation | <input type="text" value="Point to Point"/> | | | | |
| | Conc | WORK | MASTER | Lot No | |
| S1 | <input type="text" value="0"/> | <input type="text" value="8"/> | | | K <input type="text" value="N/A"/> |
| S2 | <input type="text" value="25"/> | <input type="text" value="838"/> | | | |
| S3 | <input type="text" value="50"/> | <input type="text" value="3687"/> | | | |
| S4 | <input type="text" value="87.5"/> | <input type="text" value="7471"/> | | | |
| S5 | <input type="text" value="200"/> | <input type="text" value="16027"/> | | | |
| S6 | | | | | |

Calibration (2/2) Autodilution

| | | | | |
|-----------------|--------------------------------------------|----------------------------------|---------------------------------|---------------------------------------------|
| Serial Dilution | <input checked="" type="checkbox"/> enable | <input type="checkbox"/> disable | Claculation | <input type="text" value="Point to Point"/> |
| | Conc | Post Sampling (μL) | Pre Sampling (μL) | Diluent(μL) |
| S1 | <input type="text" value="0"/> | <input type="text" value="6"/> | <input type="text" value="0"/> | <input type="text" value="80"/> |
| S2 | <input type="text" value="25"/> | <input type="text" value="6"/> | <input type="text" value="10"/> | <input type="text" value="70"/> |
| S3 | <input type="text" value="50"/> | <input type="text" value="6"/> | <input type="text" value="20"/> | <input type="text" value="60"/> |
| S4 | <input type="text" value="87.5"/> | <input type="text" value="6"/> | <input type="text" value="35"/> | <input type="text" value="45"/> |
| S5 | <input type="text" value="200"/> | <input type="text" value="6"/> | <input type="text" value="0"/> | <input type="text" value="0"/> |
| S6 | | | | |

Reagent Registration

| | | | |
|--------------|--------------------------------------------|-----------------------------------------------|-----------------------------------------------|
| Reagent Code | <input type="text" value="0173"/> | Reagent Name | <input type="text" value="RF"/> |
| R1 | <input checked="" type="checkbox"/> enable | Volume (L) <input type="text" value="**"/> mL | Volume (S) <input type="text" value="**"/> mL |
| R2 | <input checked="" type="checkbox"/> enable | Volume (L) <input type="text" value="**"/> mL | Volume (S) <input type="text" value="**"/> mL |
| | | Stability Check | Term |
| | | <input checked="" type="checkbox"/> enable | <input type="text" value="**"/> days |
| | | <input checked="" type="checkbox"/> enable | <input type="text" value="**"/> days |

** Operator definable N/A not applicable to this test Calibration curve is only as example