

*quantex* TRANSFERRIN**Kit Configuration**

P/N 3000-2172	2 x 100 mL TRF R1
	2 x 4 mL TRF R2

**Reagent Preparation**

P/N 3000-2172	TRF R1: Ready to use.
	TRF R2: Ready to use.
	Place the bottles into reagent tray.

**In Use Stability**

For optimal stability remove reagents from the system immediately after use, and store them at 2-8°C in the original vial securely closed.

**Specimen**

Serum.

**Calibration**

Use *quantex* PROTEINS standard multipoint Cat. No 300-2128. The calibrator concentrations are indicated insert sheet. Recalibrate every 60 days or when a new lot of reagent is used.

**Quality Control**

Use *quantex* PROTEINS Control I/II Cat. No. 3000-2122.

**Calculation of Analytical Results**

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the ILab 600/650 Operator's Manual.

**Reference Interval**

The reported expected range for transferrin in adults is 200-360 mg/dL (2.0-3.6 g/L).

**References / Literatur / Bibliografía / Bibliographie / Bibliografia /**

See package insert enclosed in the kit

<b>Performance Characteristics</b>
------------------------------------

**Limitation/Interfering Substances**

No significant interference from bilirubin up to concentrations of 20 mg/dL (342 µmol/L), hemoglobin up to concentrations of 500 mg/dL (0.3 mmol/L), lipemia up to concentrations of 1000 mg/dL (11.3 mmol/L) and rheumatoid factor up to 300 IU/mL. For a comprehensive review of interfering substances, refer to the publication by Young *et al.*<sup>1</sup>

**Linearity**

no rerun 40- 660 mg/dL ; with rerun 40- 6600 mg/dL

**Quantification Limit**

40 mg/dL



### Instrument Settings

<b>Photometric Test Parameters</b>		<b>Serum</b>
Test No.		**
Test Name, Test Code		TRF, TRF
Sample Type		Serum
Reporting Unit, Decimal Points		mg/dL, 1
Reaction Cycle		Standard
Twin Analysis		OFF
Methodology Type, Measuring Point		End Point, 17/33
Photometric Methodology		1 Wavelength
Primary/Secondary Wavelength		510
<b>Sampling Conditions</b>		
<i>Sampling 1</i>	Sample Vol.	3
	Sample/Diluent Vol.	0/0
<i>Sampling 2</i>	Sample Vol.	3
	Sample/Diluent Vol.	15/135
<i>Sampling 3</i>	Sample Vol.	9
	Sample/Diluent Vol.	0/0
<i>Sampling 4</i>		***
Diluent Code		Saline
Diluent Warning Limit		N/A
First Run		Sampling 1
Below/Above Normal Range		***
Panic L		***
Panic H		Sampling 2
Noise		***
Prozone		N/A
High!, ABS!		Sampling 2
Sample Volume Reduction		**
<b>Reagent Volumes</b>		
R1	Code	01811
	Rgt/Dil. Vol. Stirring	200/0, ON
	Low Vol. Warning Limit	***
	Stability (days)	**
R2	Code	01812
	Rgt/Dil. Vol. Stirring	20/20, ON
	Low Vol. Warning Limit	***
	Stability (days)	**

<b>Ranges and Evaluation Criteria</b>		<b>Serum</b>
Normal Range-Male		**
Normal Range-Female		**
Normal Range-Other		**
Valid Range		0 - 660
Hemolysis/Icterus/Lipemia Limit		***
Reaction Slope		Positive
Absorbance Limit		Above, 3200
Prozone Limit		N/A
Non Linear Limit		N/A
Slope/Intercept Correction		1/0
Qualitative Report		OFF
<b>Calibration Conditions</b>		
Calibration		5 Points, point to point, 2 Reps
Stability (days)		60
Calibrator, Concentration		Proteins Std, *
R-Blank Limit (mAbs)		500
Cal. Repts Range (%)		***
Min Cal. Response (mAbs)		***
Cal. Factor Change (%)		***
M-Point Curve Fit (%)		N/A
Reagent Blank		ON
Auto R-Blank by Bottle		ON

- \* Lot dependent
- \*\* operator definable
- \*\*\* optional
- N/A not applicable to this test