

quantex HAPTOGLOBIN

Kit Configuration

P/N 3000-2305	1 x 35 mL HAPT R1
	2 x 3.5 mL HAPT R2

Reagent Preparation

P/N 3000-2305 HAPT R1: Ready to use.
 HAPT R2: Ready to use.
 Place the bottles into reagent tray.

In Use Stability

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial securely closed.

Specimen

Serum.

Calibration

Use quantex PROTEINS standard multipoint Cat. No 3000-2128. See insert sheet for lot specific concentrations. Recalibrate every 90 days or when a new lot of reagent is used.

Quality Control

Use quantex PROTEINS Control I/II Cat. No.3000-2122.

Calculation of Analytical Results

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the ILab 600/650 Operator's Manual.

Reference Interval

The reported expected range for haptoglobin in adults is 30-200 mg/dL (0.3 – 2.0 g/L).

References / Literatur / Bibliografia / Bibliographie / Bibliografía /

See package insert enclosed in the kit

Performance Characteristics

Limitation/Interfering Substances

Non significant interference from triglycerides up to concentration of 2000 mg/dL (22.6 mmol/L). No significant interference from bilirubin up to 20 mg/dL (340 µmol/L) and hemoglobin up to concentrations of 500 mg/dL (0.3 mmol/L). For a comprehensive review of interfering substances, refer to the publication by Young *et al.*¹

Precision

	Samples/Runs	Mean (mg/dL)	CV(%)	Mean(mg/dL)	CV(%)
Within run	3/10	64.2	1.9	197.0	3.6
Total	3/10	64.2	3.7	197.0	3.8

Linearity

no rerun 5 - 260 mg/dL ; with rerun 5 - 2600 mg/dL

Instrument Settings

Photometric Test Parameters		Serum
Test No.		**
Test Name, Test Code		Hapt, Hapt
Sample Type		Serum
Reporting Unit, Decimal Points		mg/dL, 1
Reaction Cycle		Standard
Twin Analysis		OFF
Methodology Type, Measuring Point		End Point, 17/33
Photometric Methodology		1 Wavelength
Primary/Secondary Wavelength		600
Sampling Conditions		
Sampling 1	Sample Vol.	2
	Sample/Diluent Vol.	0/0
Sampling 2	Sample Vol.	2
	Sample/Diluent Vol.	15/135
Sampling 3	Sample Vol.	***
	Sample/Diluent Vol.	0/0
Sampling 4		***
Diluent Code		Saline
Diluent Warning Limit		N/A
First Run		Sampling 1
Below/Above Normal Range		***
Panic L		***
Panic H		Sampling 2
Noise		***
Prozone		N/A
High!, ABS!		Sampling 2
Sample Volume Reduction		**
Reagent Volumes		
R1	Code	01801
	Rgt/Dil. Vol. Stirring	180/0, ON
	Low Vol. Warning Limit	***
	Stability (days)	**
R2	Code	01802
	Rgt/Dil. Vol. Stirring	45/10, ON
	Low Vol. Warning Limit	***
	Stability (days)	**

Ranges and Evaluation Criteria	Serum
Normal Range-Male	**
Normal Range-Female	**
Normal Range-Other	**
Valid Range	0 - 260
Hemolysis/Icterus/Lipemia Limit	***
Reaction Slope	Positive
Absorbance Limit	Above, 3200
Prozone Limit	N/A
Non Linear Limit	N/A
Slope/Intercept Correction	1/0
Qualitative Report	OFF
Calibration Conditions	
Calibration	5 Points, point to point, 2 Reps
Stability (days)	90
Calibrator, Concentration	Protein std, *
R-Blank Limit (mAbs)	400
Cal. Reps Range (%)	***
Min Cal. Response (mAbs)	***
Cal. Factor Change (%)	***
M-Point Curve Fit (%)	N/A
Reagent Blank	ON
Auto R-Blank by Bottle	ON

* Lot dependent
 ** operator definable
 *** optional
 N/A not applicable to this test