

Kit Configuration

| | |
|---------------|--------------------------|
| P/N 3000-2307 | 1 x 40 mL β_2 m R1 |
| | 4 x 3 mL β_2 m R2 |

Reagent Preparation

| | |
|---------------|--------------------------------------|
| P/N 3000-2307 | β_2 m R1: Ready to use |
| | β_2 m R2: Ready to use. |
| | Place the bottles into reagent tray. |

In Use Stability

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial securely closed.
On board: 28 days

Specimen

Serum.

Calibration

Use quantex β_2 -microglobulin standard Cat. No 3000-2192. See vial label for lot specific concentration. A reagent blank should be run daily before sample analysis. Recalibrate every 30 days or when a new lot of reagent is used.

Quality Control

Use quantex Proteins Control I/II Cat. No. 3000-2122.

Calculation of Analytical Results

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the ILab 600/650 Operator's Manual.

Reference Interval

Serum β_2 m concentrations have been reported to be independent of body mass and sex but to be slightly increased in elderly persons. In the serum or plasma of apparently healthy persons the normal range is approximately 0.97 to 2.64 mg/L by a RIA technique. Concentrations of β_2 m in urine from healthy subjects averaged 0.098 mg/L with an upper normal limit of 0.32 mg/L.

References / Literatur / Bibliografía / Bibliographie / Bibliografia /

See package insert enclosed in the kit

Performance Characteristics

Limitation/Interfering Substances

No significant interference from lipemia up to a sample absorbance of 4.3/cm at 660 nm, bilirubin up to concentrations of 25 mg/dL (425 μ mol/L), hemoglobin up to concentrations of 400 mg/dL (0.24 mmol/L) and rheumatoid factor up to 750 IU/mL. For a comprehensive review of interfering substances, refer to the publication by Young *et al.*¹

Precision

| Serum | Samples/Runs | Mean (mg/L) | CV(%) | Mean (mg/L) | CV(%) |
|------------|--------------|-------------|-------|-------------|-------|
| Within run | 4/10 | 0.77 | 2.0 | 6.02 | 2.7 |
| Total | 4/10 | 0.77 | 2.9 | 6.02 | 3.7 |

Method Comparison

| | Serum | Urine |
|---------------------------|---------------|---------------|
| Comparison method | same reagent | same reagent |
| Comparison instrument (x) | ILab 900/1800 | ILab 900/1800 |
| Slope | 0.959 | 0.817 |
| y intercept | 0.053 | 0.003 |
| Range (mg/L) | 0.20 - 109.85 | 0.01 - 67.50 |
| Mean X (mg/L) | 7.75 | 10.02 |
| Mean Y (mg/L) | 7.48 | 8.19 |
| r | 0.997 | 0.992 |
| Syx | 1.04 | 1.75 |
| n | 62 | 40 |

Linearity

no rerun 0.5 - 16 mg/L ; with rerun 0.05 - 96 mg/L

Note:

High sensitivity application for urine samples:

Very low concentrations of β_2 m can be quantified in urine samples using the rerun capability of the instrument, when possible. Another option is using the "high sensitivity application" but only with those urine samples reporting LESS THAN 1 mg/l with the regular application, otherwise false low reported results due to prozone effect may be obtained.

Instrument Settings

| Photometric Test Parameters | | Serum | Urine High Sensitivity |
|------------------------------------|------------------------|--------------------|-------------------------------|
| Test No. | | ** | ** |
| Test Name, Test Code | | B2M, B2M | B2MHS, B2MHS |
| Sample Type | | Serum/Plasma/Urine | Urine <1 mg/L |
| Reporting Unit, Decimal Points | | mg/L, 2 | mg/L, 2 |
| Reaction Cycle | | Standard | Standard |
| Twin Analysis | | OFF | OFF |
| Methodology Type, Measuring Point | | End Point, 20/33 | End Point, 20/33 |
| Photometric Methodology | | 1 Wavelength | 1 Wavelength |
| Primary/Secondary Wavelength | | 570 | 570 |
| Sampling Conditions | | | |
| <i>Sampling 1</i> | Sample Vol. | 3 | 30 |
| | Sample/Diluent Vol. | 0/0 | 0/0 |
| <i>Sampling 2</i> | Sample Vol. | 3 | *** |
| | Sample/Diluent Vol. | 30/150 | *** |
| <i>Sampling 3</i> | Sample Vol. | 30 | *** |
| | Sample/Diluent Vol. | 0/0 | *** |
| <i>Sampling 4</i> | | *** | *** |
| Diluent Code | | Saline | Saline |
| Diluent Warning Limit | | *** | *** |
| First Run | | Sampling 1 | Sampling 1 |
| Below/Above Normal Range | | *** | *** |
| Panic L | | Sampling 3 | *** |
| Panic H | | Sampling 2 | *** |
| Noise | | *** | *** |
| Prozone | | N/A | N/A |
| High!, ABS! | | Sampling 2 | *** |
| Sample Volume Reduction | | ** | ** |
| Reagent Volumes | | | |
| R1 | Code | 01821 | 01821 |
| | Rgt/Dil. Vol. Stirring | 150/0, ON | 150/0, ON |
| | Low Vol. Warning Limit | *** | *** |
| | Stability (days) | 28 | ** |
| R2 | Code | 01822 | 01822 |
| | Rgt/Dil. Vol. Stirring | 120/20, ON | 120/20, ON |
| | Low Vol. Warning Limit | *** | *** |
| | Stability (days) | 28 | ** |

| Ranges and Evaluation Criteria | Serum | Urine High Sensitivity |
|---------------------------------------|-------------------------|-------------------------------|
| Normal Range-Male | ** | ** |
| Normal Range-Female | ** | ** |
| Normal Range-Other | ** | ** |
| Valid Range | 0.50 - 16.00 | 0.05 - 1.60 |
| Hemolysis/Icterus/Lipemia Limit | *** | *** |
| Reaction Slope | Positive | Positive |
| Absorbance Limit | Above, 3200 | Above, 3200 |
| Prozone Limit | N/A | N/A |
| Non Linear Limit | N/A | N/A |
| Slope/Intercept Correction | 1/0 | 1/0 |
| Qualitative Report | OFF | OFF |
| Calibration Conditions | | |
| Calibration | 1 Point, Linear, 3 Reps | 1 Point, Linear, 3 Reps |
| Stability (days) | 30 | ** |
| Calibrator, Concentration | Std β2m * | Std β2m * |
| R-Blank Limit (mAbs) | 1000 | 1000 |
| Cal. Reps Range (%) | *** | *** |
| Min Cal. Response (mAbs) | *** | *** |
| Cal. Factor Change (%) | *** | *** |
| M-Point Curve Fit (%) | N/A | N/A |
| Reagent Blank | ON | ON |
| Auto R-Blank by Bottle | ON | ON |

* Lot dependent
 ** operator definable
 *** optional
 N/A not applicable to this test

Note: for urine High Sensitivity application the Calibrator has to be dilute 1:10 (1+9) with Saline. The Calibrator concentration is calculated dividing Calibrator value by 10.