

quantex HbA_{1c}

Kit Configuration

P/N 3000-2314	1 x 50 mL HbA _{1c} (Diluent)
	1 x 28 mL Hb R1
	1 x 14 mL HbA _{1c} R1
	1 x 14 mL HbA _{1c} R2

Reagent Preparation

P/N 3000-2314	HbA _{1c} (Diluent):	Ready to use
	Hb R1:	Ready to use. Invert to mix well before first use. Avoid foam formation.
	HbA _{1c} R1:	Ready to use. Invert to mix well before first use. Avoid foam formation.
	mL HbA _{1c} R2:	Ready to use. Invert to mix well before first use. Avoid foam formation.

Place the bottles into reagent tray.

In Use Stability

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial securely closed.

Specimen

Whole blood or capillary blood samples.

Calibration

Use quantex HbA_{1c} standard multipoint Cat. No 3000-2315.

Calibration of the quantex Hb hemoglobin reagent is performed with quantex HbA_{1c} standard level 1 only.

Calibration of the quantex HbA_{1c} reagents R1 and R2 is performed with quantex HbA_{1c} standards levels 1 to 6.

The concentrations are indicated in the insert sheet. A reagent blank should be run daily before sample analysis. Recalibrate every 20 days or when a new lot of reagent is used.

Quality Control

Use quantex HbA_{1c} Control I/II Cat. No.3000-2316.

Calculation of Analytical Results

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the ILab 600/650 Operator's Manual.

Reference Interval

Depending on the assay used, HbA_{1c} is approximately 4-6% in non diabetics, 6-8% in controlled diabetics, and can be as much as 20% in uncontrolled diabetics. 124 apparently normal healthy donors undergoing physical examination were tested for HbA_{1c} using the quantex HbA_{1c} reagents. The range of results of HbA_{1c} was of 4.5-6.2% and a mean of 5.4% resulted.

In any case, it is recommended that each laboratory established its own expected range.

References / Literatur / Bibliografia / Bibliographie / Bibliografia /

See package insert enclosed in the kit

Performance Characteristics

Limitation/Interfering Substances

No significant interference from triglycerides up to concentrations of 1600 mg/dL (18 mmol/L) , bilirubin up to concentrations of 30 mg/dL (513 μmol/L) and rheumatoid factor up to 2000 IU/mL. For a comprehensive review of interfering substances, refer to the publication by Young *et al.*¹

Linearity

Total Hb : 7 to 23 g/dL

% HbA_{1c} : 1.8% to 17.2% (for total Hb of 14 g/dL)

Samples with values above 17.2% of HbA_{1c} should not be diluted and results should be reported as > 17.2% HbA_{1c}.

Instrument Settings

Photometric Test Parameters		
	HbA_{1c}	Hb Total
Test No.	**	**
Test Name, Test Code	HbA _{1c} ,HbA _{1c}	THb
Sample Type	Others	Others
Reporting Unit, Decimal Points	g/dL, 2	g/dL, 2
Reaction Cycle	Standard	Standard
Twin Analysis	OFF	OFF
Methodology Type, Measuring Point	End Point, 19/26	End Point, 17
Photometric Methodology	1 Wavelength	1 Wavelength
Primary/Secondary Wavelength	700	600
Sampling Conditions		
Sampling 1	Sample Vol. 4	15
	Sample/Diluent Vol. 0/0	0/0
Sampling 2	Sample Vol. 2	5
	Sample/Diluent Vol. 0/0	0/0
Sampling 3	Sample Vol. 8	30
	Sample/Diluent Vol. 0/0	0/0
Sampling 4	***	***
Diluent Code	Saline	Saline
Diluent Warning Limit	***	***
First Run	Sampling 1	Sampling 1
Below/Above Normal Range	***	***
Panic L	***	***
Panic H	Sampling 2	Sampling 2
Noise	***	***
Prozone	N/A	N/A
High!, ABS!	Sampling 2	Sampling 2
Sample Volume Reduction	**	**
Reagent Volumes		
R1	Code 01981	01971
	Rgt/Dil. Vol. Stirring 95/5, ON	195/5, ON
	Low Vol. Warning Limit ***	***
	Stability (days) **	**
R2	Code 01982	N/A
	Rgt/Dil. Vol. Stirring 95/5, ON	N/A
	Low Vol. Warning Limit ***	N/A
	Stability (days) **	N/A

Ranges and Evaluation Criteria		
	HbA_{1c}	Hb Total
Normal Range-Male	**	**
Normal Range-Female	**	**
Normal Range-Other	**	**
Valid Range	0 - 2.4	0-23
Hemolysis/Icterus/Lipemia Limit	***	***
Reaction Slope	Positive	Positive
Absorbance Limit	Above, 3200	Above, 3200
Prozone Limit	N/A	N/A
Non Linear Limit	N/A	N/A
Slope/Intercept Correction	1/0	1/0
Qualitative Report	OFF	OFF
Calibration Conditions		
Calibration	6 Points, point to point, 2 Reps	1 Point, Linear, 3 Reps
Stability (days)	20	***
Calibrator, Concentration	HbA _{1c} Std, *	Hb Std, *
R-Blank Limit (mAbs)	1800	100
Cal. Reps Range (%)	***	***
Min Cal. Response (mAbs)	***	***
Cal. Factor Change (%)	***	***
M-Point Curve Fit (%)	N/A	N/A
Reagent Blank	ON	ON
Auto R-Blank by Bottle	ON	ON

Calculated Test Parameters	
Test n°	**
Test Name, Test Code	%HbA _{1c} , %HbA _{1c}
Sample Type	N/A
Reporting Unit, Decimal Points	%, 2
Equation	$([HbA_{1c}] / [Hb]) * 100$
Normal Range Male	**
Normal Range Female	**
Normal Range Others	**
Valid Range	1.8 - 17.2

* Lot dependent
 ** operator definable
 *** optional
 N/A not applicable to this test