

## quantex CARBAMAZEPINE

**Kit Configuration**

P/N 3000-2277	2 x 11 mL CARB R1
	2 x 3 mL CARB R2

**Reagent Preparation**

P/N 3000-2277: CARB R1 (Buffer) : Ready to use.  
 CARB R2 (Reagent) : Ready to use. Invert to mix well before first use. Avoid foam formation  
 Place the bottles into reagent tray.

**In Use Stability**

Until the expiration date shown on the vial when stored at 2-8°C. For optimal stability remove reagents from the system and store them at 2-8°C in the original vial securely closed.

**Specimen**

Serum,.

**Calibration**

Use quantex CARBAMAZEPINE standard multipoint Cat. No 3000-2318. The concentrations in µg/mL are indicated on the vial labels. Recalibrate every 90 days, when a new lot of reagents is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

**Quality Control**

Use quantex TDM control I/II Cat. No 3000-2303.

**Calculation of Analytical Results**

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the ILab 350 Operator Manual.

**Therapeutic Range**

The typical therapeutic range is 4 - 10 µg/mL (17 - 42 µmol/L) and the toxic range is > 15 µg/mL (> 63 µmol/L).

To convert results to µmol/L multiply by 4.23.

**References / Literatur / Bibliografía / Bibliographie / Bibliografia /**

See package insert enclosed in the kit

**Performance Characteristics****Limitation/Interfering Substances**

No significant interference from bilirubin up to concentrations of 20 mg/dL, hemoglobin up to concentrations of 1000 mg/dL and lipemia up to concentrations of 20 g/L. For a comprehensive review of interfering substances, refer to the publication by Young *et al.*<sup>1</sup>

**Precision**

	Samples/ Runs	Mean (µg/mL)	CV (%)	Mean (µg/mL)	CV (%)
Within run	4/10	3.2	3.7	8.5	7.7
Total	4/10	3.2	2.9	8.5	7.7

**Linearity**

no rerun 1.36 to 20 µg/mL

With rerun 1.36 to 80 µg/mL



### Instrument Settings

<b>Chemistry Parameters</b>				<b>R1</b>			
Method	<input type="text"/>	Reagent Name	<input type="text" value="CBZ"/>	Volume	<input type="text" value="220 μL"/>		
Name	<input type="text" value="CBZ"/>	R2	<input type="text" value="enable"/>				
Unit	<input type="text" value="μg/mL"/>	Reagent Name	<input type="text" value="CBZ"/>	Volume	<input type="text" value="60 μL"/>		
Assay Type	<input type="text" value="End"/>	Wash	<input type="text" value="disable"/>	Reagent Name	<input type="text"/>		
				Diluent	<input type="text" value="enable"/>	Reagent Type	<input type="text" value="Saline"/>
<b>Measuring Points</b>				<b>Decimal Points</b>			
	<input type="text" value="1"/>	<input type="text" value="enable"/>	start	<input type="text" value="14"/>			
			end	<input type="text" value="15"/>			
	<input type="text" value="2"/>	<input type="text" value="enable"/>	start	<input type="text" value="25"/>	<b>Normal Range</b>		
			end	<input type="text" value="26"/>	<input type="text" value="4"/>	<input type="text" value="10"/>	
<b>Wave Length</b>				<b>Technical Range (Conc)</b>			
Prim	<input type="text" value="600"/>	Sec	<input type="text"/>	<input type="text" value="0.0"/>		<input type="text" value="20"/>	
				mAbs/10 <input type="text" value="-30000 / 30000"/>			
<b>Sampling Volume</b>				<b>RPT Wash</b>			
Dilution	<input type="text" value="2 μL"/>	<input type="text" value="disable"/>		(R1)	<input type="text" value="Sys Water"/>		
	<input type="text"/>	<input type="text"/>	<input type="text"/>	(R2)	<input type="text" value="Sys Water"/>		
<b>Rerun ( High )</b>				<b>Instrument Factor a</b>			
Dilution	<input type="text" value="2 μL"/>	<input type="text" value="enable"/>		<input type="text" value="1"/>	<b>b</b>		<input type="text" value="0"/>
	<input type="text" value="35 μL"/>	<input type="text" value="105 μL"/>	<b>Stirring Speed</b>		R1	<input type="text" value="high"/>	R2
<b>Rerun ( Low )</b>							
	<input type="text" value="4 μL"/>						

### Calibration Checks

<b>** Duplicate Limit</b>	<input type="text" value="**"/>	mAbs/10	<b>Sampling Method for Standards</b>				
<b>** Sensitivity Limit</b>	<input type="text" value="**"/>	mAbs/10	<input checked="" type="checkbox"/>	Duplicate			
			<input type="checkbox"/>	Triplicate			
<b>** Linearity Limit</b>	<input type="text" value="**"/>	%	<b>Blank measurement</b>				
<b>** Prozone Limit</b>	<input type="text" value="**"/>	upper	<input checked="" type="checkbox"/>	Enable Reagent blank			
SL1-S	<input type="text" value="**"/>	SL1-F	<input type="text" value="None"/>				
SL2-S	<input type="text" value="**"/>	SL2-F	<b>Reagent blank measurement at calibration</b>				
Sens	<input type="text" value="**"/>	mAbs/10	<input checked="" type="checkbox"/>	Reagent blank (system water)			
<input checked="" type="checkbox"/>	<b>Absorbance Limit</b>		<input type="text" value="**"/>	<b>Multiplex measurement is the same as standards</b>			
Reaction	<input type="text" value="Increase"/>		<b>Reagent Blank Limit Checks</b>				
Limit	<input type="text" value="25000"/>		<input type="text" value="mAbs/10"/>	<input type="text" value="**"/>	Duplicate limit		
			<input type="text" value="50"/>	mAbs/10			

### Calibration

<b>Method</b>	<input type="text"/>	<b>Name</b>	<input type="text" value="CBZ"/>	<b>Interval</b>	<input type="text" value="90"/>	days
<b>Calculation</b>	<input type="text" value="Point to Point"/>					
	Conc	WORK	MASTER	Lot No		
S1	<input type="text" value="0"/>	<input type="text" value="-172"/>	<input type="text"/>	<input type="text"/>	K	<input type="text" value="N/A"/>
S2	<input type="text" value="2"/>	<input type="text" value="-2063"/>	<input type="text"/>	<input type="text"/>		
S3	<input type="text" value="4"/>	<input type="text" value="-3339"/>	<input type="text"/>	<input type="text"/>		
S4	<input type="text" value="8"/>	<input type="text" value="-4661"/>	<input type="text"/>	<input type="text"/>		
S5	<input type="text" value="12"/>	<input type="text" value="-5374"/>	<input type="text"/>	<input type="text"/>		
S6	<input type="text" value="20"/>	<input type="text" value="-5214"/>	<input type="text"/>	<input type="text"/>		

### Reagent Registration

Reagent Code	<input type="text" value="0165"/>													
Reagent Name	<input type="text" value="CBZ"/>													
R1	<input checked="" type="checkbox"/>	enable	Volume (L)	<input type="text" value="**"/>	mL	Volume (S)	<input type="text" value="**"/>	mL	Stability Check	<input checked="" type="checkbox"/>	enable	Term	<input type="text" value="**"/>	days
R2	<input checked="" type="checkbox"/>	enable	Volume (L)	<input type="text" value="**"/>	mL	Volume (S)	<input type="text" value="**"/>	mL	Stability Check	<input checked="" type="checkbox"/>	enable	Term	<input type="text" value="**"/>	days
** Operator definable			N/A not applicable to this test			Calibration curve is only as example								