

quantex Valproic Acid

**Kit Configuration**

P/N 3000-2280	2 x 15 mL VALP R1
	2 x 3.8 mL VALP R2

**Reagent Preparation**

P/N 3000-2280: VALP R1: Ready to use.  
 VALP R2: Ready to use. Invert to mix well before first use. Avoid foam formation  
 Place the bottles into reagent tray.

**In Use Stability**

Stable until the expiration date shown on the vial when stored at 2-8°C. For optimal stability remove reagents from the system and store them at 2-8°C in the original vial securely closed.

**Specimen**

Serum,.

**Calibration**

Use quantex VALPROIC ACID standard multipoint Cat. No 3000-2288. The concentrations in µg/mL are indicated on the vial labels. Recalibrate every 60 days, when a new lot of reagents is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

**Quality Control**

Use quantex TDM control I/II Cat. No 3000-2303.

**Calculation of Analytical Results**

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the I Lab 350 Operator Manual.

**Therapeutic Range**

The typical therapeutic range is 50 - 100 µg/mL (347 - 693 µmol/L) and the toxic range is > 100 µg/mL (> 693 µmol/L). Some patients achieve the desired therapeutic response at levels outside this range; therefore, individual clinical evaluation should be considered when interpreting assay results.  
 To convert results to µmol/l multiply by 6.93.

**References / Literatur / Bibliografía / Bibliographie / Bibliografia /**

See package insert enclosed in the kit

**Performance Characteristics**

**Limitation/Interfering Substances**

No significant interference from bilirubin up to concentrations of 20 mg/dL, hemoglobin up to concentrations of 1000 mg/dL and lipemia up to concentrations of 20 g/L. For a comprehensive review of interfering substances, refer to the publication by Young *et al.*<sup>1</sup>

**Precision**

	Samples/ Runs	Mean (µg/mL)	CV (%)	Mean (µg/mL)	CV (%)
Within run	4/10	37.8	2.7	88.8	1.1
Total	4/10	37.8	6.4	88.8	5.1

**Linearity**

no rerun 11 to 150 µg/mL  
 With rerun 11 to 600 µg/mL

### Instrument Settings

<b>Chemistry Parameters</b>				<b>R1</b>			
Method	<input type="text"/>	Reagent Name	<input type="text" value="VLP"/>	Volume	<input type="text" value="200 μL"/>		
Name	<input type="text" value="VLP"/>	R2	<input type="text" value="enable"/>				
Unit	<input type="text" value="μg/mL"/>	Reagent Name	<input type="text" value="VLP"/>	Volume	<input type="text" value="50 μL"/>		
Assay Type	<input type="text" value="End"/>	Wash	<input type="text" value="disable"/>	Reagent Name	<input type="text"/>		
				Diluent	<input type="text" value="enable"/>	Reagent Type	<input type="text" value="Saline"/>
<b>Measuring Points</b>				<b>Decimal Points</b>			
	<input type="text" value="1"/>	<input type="text" value="enable"/>	start	<input type="text" value="14"/>			
			end	<input type="text" value="15"/>			
	<input type="text" value="2"/>	<input type="text" value="enable"/>	start	<input type="text" value="25"/>	<b>Normal Range</b>		
			end	<input type="text" value="26"/>	<input type="text" value="50"/>	<input type="text" value="100"/>	
<b>Wave Length</b>				<b>Technical Range (Conc)</b>			
Prim	<input type="text" value="600"/>	Sec	<input type="text"/>	<input type="text" value="0.0"/>		<input type="text" value="150"/>	
				mAbs/10 <input type="text" value="-30000 / 30000"/>			
<b>Sampling Volume</b>				<b>RPT Wash</b>			
Dilution	<input type="text" value="2 μL"/>	<input type="text" value="disable"/>		(R1)	<input type="text" value="Sys Water"/>		
	<input type="text"/>	<input type="text"/>		(R2)	<input type="text" value="Sys Water"/>		
<b>Rerun ( High )</b>				<b>Instrument Factor a</b>			
Dilution	<input type="text" value="2 μL"/>	<input type="text" value="enable"/>		<input type="text" value="1"/>	<b>b</b>		<input type="text" value="0"/>
	<input type="text" value="35 μL"/>	<input type="text" value="105 μL"/>		<b>Stirring Speed</b>			
<b>Rerun ( Low )</b>				R1	<input type="text" value="high"/>	R2	<input type="text" value="high"/>

### Calibration Checks

**	<b>Duplicate Limit</b>	<input type="text"/>	**	mAbs/10	<b>Sampling Method for Standards</b>		
**	<b>Sensitivity Limit</b>	<input type="text"/>	**	mAbs/10	<input checked="" type="checkbox"/>	Duplicate	
					<input type="checkbox"/>	Triplicate	
**	<b>Linearity Limit</b>	<input type="text"/>	**	%	<b>Blank measurement</b>		
**	<b>Prozone Limit</b>	<input type="text"/>	upper		<input checked="" type="checkbox"/>	Enable Reagent blank	
	SL1-S	**	SL1-F	**	<input type="text" value="None"/>		
	SL2-S	**	SL2-F	**	<b>Reagent blank measurement at calibration</b>		
	Sens	<input type="text"/>	**	mAbs/10	<input checked="" type="checkbox"/>	Reagent blank (system water)	
<input checked="" type="checkbox"/>	<b>Absorbance Limit</b>	<input type="text"/>	**	Multiplex measurement is the same as standards	<b>Reagent Blank Limit Checks</b>		
	Reaction Limit	<input type="text" value="Increase"/>	<input type="text" value="25000"/>	mAbs/10	**	Duplicate limit <input type="text" value="50"/>	

### Calibration

<b>Method</b>	<input type="text"/>	<b>Name</b>	<input type="text" value="VLP"/>	<b>Interval</b>	<input type="text" value="60"/>	days
<b>Calculation</b>	<input type="text" value="Point to Point"/>					
	Conc	WORK	MASTER	Lot No		
S1	<input type="text" value="0"/>	<input type="text" value="-632"/>	<input type="text"/>	<input type="text"/>	K	<input type="text" value="N/A"/>
S2	<input type="text" value="12.5"/>	<input type="text" value="-1833"/>	<input type="text"/>	<input type="text"/>		
S3	<input type="text" value="25"/>	<input type="text" value="-3240"/>	<input type="text"/>	<input type="text"/>		
S4	<input type="text" value="50"/>	<input type="text" value="-5544"/>	<input type="text"/>	<input type="text"/>		
S5	<input type="text" value="100"/>	<input type="text" value="-7241"/>	<input type="text"/>	<input type="text"/>		
S6	<input type="text" value="150"/>	<input type="text" value="-7506"/>	<input type="text"/>	<input type="text"/>		

### Reagent Registration

Reagent Code	<input type="text" value="0150"/>													
Reagent Name	<input type="text" value="VLP"/>													
R1	<input checked="" type="checkbox"/>	enable	Volume (L)	<input type="text"/>	mL	Volume (S)	<input type="text"/>	mL	Stability Check	<input checked="" type="checkbox"/>	enable	Term	<input type="text"/>	days
R2	<input checked="" type="checkbox"/>	enable	Volume (L)	<input type="text"/>	mL	Volume (S)	<input type="text"/>	mL	Stability Check	<input checked="" type="checkbox"/>	enable	Term	<input type="text"/>	days
**	Operator definable		N/A		not applicable to this test			Calibration curve is only as example						