

quantex AMIKACIN



Kit Configuration

P/N 3000-2317	2 x 12.5 mL AMIK R1
	2 x 3.8 mL AMIK R2

Reagent Preparation

P/N 3000-2317: AMIK R1: Ready to use.
 AMIK R2: Ready to use. Invert to mix well before first use. Avoid foam formation
 Place the bottles into reagent tray.

In Use Stability

For optimal stability remove reagents from the system and store them at 2-8°C in the original vial securely closed.

Specimen

Serum,

Calibration

Use quantex AMIKACIN standard multipoint Cat. No 3000-2318. The concentrations in µg/mL are indicated on the vial labels. Recalibrate every 36 days, when a new lot of reagents is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

Quality Control

Use quantex TDM control I/II Cat. No 3000-2303.

Calculation of Analytical Results

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the ILab 350 Operator Manual.

Therapeutic Range

Recommended "trough" levels are: 5 - 8 µg/mL (8.5 - 14 µmol/L).

Recommended "peak" levels are: 20 - 30 µg/mL (34 - 51 µmol/L).

"Trough" concentrations above 10 µg/mL (> 17 µmol/L) and "peak" concentrations above 35 µg/mL (> 60 µmol/L) are often associated with renal impairment and ototoxicity.

To convert results to µmol/L multiply by 1.71.

References / Literatur / Bibliografia / Bibliographie / Bibliografia /

See package insert enclosed in the kit

Performance Characteristics

Limitation/Interfering Substances

No significant interference from bilirubin up to concentrations of 20 mg/dL, hemoglobin up to concentrations of 1000 mg/dL and lipemia up to concentrations of 20 g/L. For a comprehensive review of interfering substances, refer to the publication by Young *et al.*¹

Precision

	Samples/ Runs	Mean (µg/mL)	CV (%)	Mean (µg/mL)	CV (%)
Within run	4/10	5.1	7.1	16.6	2.6
Total	4/10	5.1	8.3	16.6	3.0

Linearity

no rerun 1.5 to 50 µg/mL

With rerun 1.5 to 75 µg/mL



Instrument Settings

Chemistry Parameters				R1			
Method Name	<input type="text" value="AMIK"/>	Reagent Name	<input type="text" value="AMIK"/>	Volume	<input type="text" value="250 μL"/>		
Unit	<input type="text" value="μg/mL"/>	Reagent Name	<input type="text" value="AMIK"/>	Volume	<input type="text" value="75 μL"/>		
Assay Type	<input type="text" value="End"/>	Wash	<input type="text" value="disable"/>	Reagent Name			
		Diluent	<input type="text" value="disable"/>	Reagent Type			
Measuring Points	1 enable	start	<input type="text" value="14"/>	Decimal Points	<input type="text" value="1"/>		
		end	<input type="text" value="15"/>	Normal Range	<input type="text" value="0"/>	<input type="text" value="50"/>	
	2 enable	start	<input type="text" value="25"/>				
		end	<input type="text" value="26"/>				
Wave Length Prim	<input type="text" value="700"/>	Sec	<input type="text" value=""/>	Technical Range (Conc) mAbs/10	<input type="text" value="0.0"/>	<input type="text" value="60"/>	
Sampling Volume	<input type="text" value="3 μL"/>	RPT Wash (R1)	<input type="text" value="Sys Water"/>				
Dilution	<input type="text" value="disable"/>	RPT Wash (R2)	<input type="text" value="Sys Water"/>				
Rerun (High)	<input type="text" value="2 μL"/>	Instrument Factor a	<input type="text" value="1"/>	b	<input type="text" value="0"/>		
Dilution	<input type="text" value="disable"/>	Stirring Speed R1	<input type="text" value="mid"/>	R2	<input type="text" value="mid"/>		
Rerun (Low)	<input type="text" value="6 μL"/>						

Calibration Checks

** Duplicate Limit	<input type="text" value="**"/>	mAbs/10	Sampling Method for Standards			
** Sensitivity Limit	<input type="text" value="**"/>	mAbs/10	<input checked="" type="checkbox"/> Duplicate			
** Linearity Limit	<input type="text" value="**"/>	%	<input type="checkbox"/> Triplicate			
** Prozone Limit	<input type="text" value="**"/>	upper	Blank measurement			
SL1-S	<input type="text" value="**"/>	SL1-F	<input checked="" type="checkbox"/> Enable Reagent blank			
SL2-S	<input type="text" value="**"/>	SL2-F	<input type="text" value="None"/>			
Sens	<input type="text" value="**"/>	mAbs/10	Reagent blank measurement at calibration			
<input checked="" type="checkbox"/> Absorbance Limit	<input type="text" value="Increase"/>		<input checked="" type="checkbox"/> Reagent blank (system water)			
Reaction Limit	<input type="text" value="25000"/>	mAbs/10	Multiplex measurement is the same as standards			
			Reagent Blank Limit Checks			
			<input type="text" value="**"/>	Duplicate limit	<input type="text" value="50"/>	mAbs/10

Calibration

Method	<input type="text" value=""/>	Name	<input type="text" value="AMIK"/>	Interval	<input type="text" value="36"/>	days
Calculation	<input type="text" value="Point to Point"/>					
	Conc	WORK	MASTER	Lot No		
S1	<input type="text" value="0"/>	<input type="text" value="54"/>			K	<input type="text" value="N/A"/>
S2	<input type="text" value="3"/>	<input type="text" value="-896"/>				
S3	<input type="text" value="10"/>	<input type="text" value="-2557"/>				
S4	<input type="text" value="20"/>	<input type="text" value="-4090"/>				
S5	<input type="text" value="35"/>	<input type="text" value="-5045"/>				
S6	<input type="text" value="50"/>	<input type="text" value="-5377"/>				

Reagent Registration

Reagent Code	<input type="text" value="199"/>											
Reagent Name	<input type="text" value="AMIK"/>											
R1	<input checked="" type="checkbox"/> enable	Volume (L)	<input type="text" value="**"/>	mL	Volume (S)	<input type="text" value="**"/>	mL	Stability Check	<input checked="" type="checkbox"/> enable	Term	<input type="text" value="**"/>	days
R2	<input checked="" type="checkbox"/> enable	Volume (L)	<input type="text" value="**"/>	mL	Volume (S)	<input type="text" value="**"/>	mL	Stability Check	<input checked="" type="checkbox"/> enable	Term	<input type="text" value="**"/>	days
**	Operator definable		N/A not applicable to this test		Calibration curve is only as example							