

## quantex AMIKACIN

### Kit Configuration

P/N 3000-2317	2 x 12.5 mL AMIK R1
	2 x 3.8 mL AMIK R2

### Reagent Preparation

P/N 3000-2317: AMIK R1: Ready to use.  
 AMIK R2: Ready to use. Invert to mix well before first use. Avoid foam formation  
 Place the bottles into reagent tray.

### In Use Stability

Stable until the expiration date shown on the vial when stored at 2-8°C. For optimal stability remove reagents from the system and store them at 2-8°C in the original vial securely closed.

### Specimen

Serum,.

### Calibration

Use quantex AMIKACIN standard multipoint Cat. No 3000-2318. The concentrations in µg/mL are indicated on the vial labels. Recalibrate every 7 days, when a new lot of reagents is used, when control recovery falls out of the expected range or when adjustments are made to the instrument. A reagent blank should be run daily before sample analysis.

### Quality Control

Use quantex TDM control I/II Cat. No 3000-2303.

### Calculation of Analytical Results

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the ILab 600/650 Operator's Manual.

### Therapeutic Range

Recommended "trough" levels are: 5 - 8 µg/mL (8.5 - 14 µmol/L).

Recommended "peak" levels are: 20 - 30 µg/mL (34 - 51 µmol/L).

"Trough" concentrations above 10 µg/mL (> 17 µmol/l) and "peak" concentrations above 35 µg/mL (> 60 µmol/l) are often associated with renal impairment and ototoxicity.

To convert results to µmol/L multiply by 1.71.<sup>1,5</sup>

### References / Literatur / Bibliografía / Bibliographie / Bibliografia /

See package insert enclosed in the kit

## Performance Characteristics

### Limitation/Interfering Substances

No significant interference from bilirubin up to concentrations of 20 mg/dL, hemoglobin up to concentrations of 1000 mg/dL and lipemia up to concentrations of 20 g/L. For a comprehensive review of interfering substances, refer to the publication by Young *et al.*<sup>4</sup>

### Precision

	Samples/ Runs	Mean (µg/mL)	CV (%)	Mean (µg/mL)	CV (%)
Within run	3/10	4.6	10.0	15.7	2.7
Total	3/10	4.6	9.2	15.7	3.4

### Linearity

no rerun 3.0 to 50 µg/mL

With rerun 3.0 to 400 µg/mL

## Instrument Settings

<b>Photometric Test Parameters</b>		<b>Serum</b>
Test No.		**
Test Name, Test Code		AMIK, AMIK
Sample Type		Serum
Reporting Unit, Decimal Points		µg/mL, 1
Reaction Cycle		Standard
Twin Analysis		OFF
Methodology Type, Measuring Point		End Point, 20/33
Photometric Methodology		1 Wavelength
Primary/Secondary Wavelength		700
<b>Sampling Conditions</b>		
<i>Sampling 1</i>	Sample Vol.	3
	Sample/Diluent Vol.	0/0
<i>Sampling 2</i>	Sample Vol.	3
	Sample/Diluent Vol.	20/140
<i>Sampling 3</i>	Sample Vol.	6
	Sample/Diluent Vol.	0/0
<i>Sampling 4</i>		***
Diluent Code		Saline
Diluent Warning Limit		***
First Run		Sampling 1
Below/Above Normal Range		***
Panic L		***
Panic H		Sampling 2
Noise		***
Prozone		N/A
High!, ABS!		Sampling 2
Sample Volume Reduction		**
<b>Reagent Volumes</b>		
R1	Code	01991
	Rgt/Dil. Vol. Stirring	250/0, ON
	Low Vol. Warning Limit	***
	Stability (days)	**
R2	Code	01992
	Rgt/Dil. Vol. Stirring	75/10, ON
	Low Vol. Warning Limit	***
	Stability (days)	**

<b>Ranges and Evaluation Criteria</b>	<b>Serum</b>
Normal Range-Male	**
Normal Range-Female	**
Normal Range-Other	**
Valid Range	0 - 50
Hemolysis/Icterus/Lipemia Limit	***
Reaction Slope	Positive
Absorbance Limit	Above, 3200
Prozone Limit	N/A
Non Linear Limit	N/A
Slope/Intercept Correction	1/0
Qualitative Report	OFF
<b>Calibration Conditions</b>	
Calibration	6 Points, point to point, 2 Reps
Stability (days)	7
Calibrator, Concentration	Amikacin Std *
R-Blank Limit (mAbs)	3500
Cal. Repts Range (%)	***
Min Cal. Response (mAbs)	***
Cal. Factor Change (%)	***
M-Point Curve Fit (%)	***
Reagent Blank	ON
Auto R-Blank by Bottle	ON

- \* Lot dependent
- \*\* operator definable
- \*\*\* optional
- N/A not applicable to this test

