

DRI® Primidone Assay**Kit Configuration**

P/N W150291	1 x 25 mL Antibody/substrate R1
	1 x 8 mL Enzyme Conjugate R2

Reagent Preparation

P/N 150291:	PRIM R1: Ready to use.
	PRIM R2: Ready to use.
	Place the bottles into reagent tray.

In Use Stability

Stable until the expiration date shown on the vial when stored at 2-8°C. For optimal stability remove reagents from the system and store them at 2-8°C in the original vial securely closed.

Specimen

Serum or plasma.

Calibration

Use Primidone calibrator multipoint Cat. No 15100185. Recalibrate the instrument when new reagents are used. A reagent blank should be run daily before sample analysis.

Quality Control

Use TDM control 3 levels. Cat. No 15755095.

Calculation of Analytical Results

The results concentration is automatically calculated by the instrument against the Calibration curve. For detailed description, refer to the Instrument settings and to the ILab 600/650 Operator's Manual.

Therapeutic Range

In most patients, a serum primidone therapeutic range of 5-12 µg/mL (17-51 µmol/L) can effectively suppress seizures without apparent side effects. Serum primidone concentrations above 15 µg/mL are often associated with toxicity.² The therapeutic range is provided only as a general guide.

References / Literatur / Bibliografia / Bibliographie / Bibliografía /

See package insert enclosed in the kit

Performance Characteristics**Limitation/Interfering Substances**

No interference with hemolyzed, lipemic and icteric samples was found with hemoglobin, cholesterol and bilirubin up to a concentration of 800 mg/dL, 400 mg/dL and 34 mg/dL, respectively.

For a comprehensive review of interfering substances, refer to the publication by Young *et al.*¹



Instrument Settings

Photometric Test Parameters		Serum
Test No.		**
Test Name, Test Code		PRIM, PRIM
Sample Type		Serum
Reporting Unit, Decimal Points		µg/mL, 2
Reaction Cycle		Standard
Twin Analysis		OFF
Methodology Type, Measuring Point		Rate, 21/25
Photometric Methodology		2 Wavelength
Primary/Secondary Wavelength		340 - 510
Sampling Conditions		
<i>Sampling 1</i>	Sample Vol.	10
	Sample/Diluent Vol.	0/0
<i>Sampling 2</i>	Sample Vol.	5
	Sample/Diluent Vol.	0/0
<i>Sampling 3</i>	Sample Vol.	20
	Sample/Diluent Vol.	0/0
<i>Sampling 4</i>		0
Diluent Code		Water
Diluent Warning Limit		***
First Run		Sampling 1
Below/Above Normal Range		***
Panic L		Sampling 3
Panic H		Sampling 2
Noise		Sampling 1
Prozone		N/A
High!, ABS!		Sampling 2
Sample Volume Reduction		**
Reagent Volumes		
R1	Code	PRI 1
	Rgt/Dil. Vol. Stirring	210/0, ON
	Low Vol. Warning Limit	0
	Stability (days)	**
R2	Code	PRI 2
	Rgt/Dil. Vol. Stirring	70/0, ON
	Low Vol. Warning Limit	0
	Stability (days)	**

Ranges and Evaluation Criteria	Serum
Normal Range-Male	5 - 12
Normal Range-Female	5 - 12
Normal Range-Other	5 - 12
Valid Range	1 - 24
Hemolysis/Icterus/Lipemia Limit	***
Reaction Slope	Positive
Absorbance Limit	Above, 3200
Prozone Limit	Above, 0, None, 0
Non Linear Limit	0
Slope/Intercept Correction	1/0
Qualitative Report	OFF
Calibration Conditions	
Calibration	Multi Points,4PLog,2 Reps
Stability (days)	**
Calibrator, Concentration	Primi Std, *
R-Blank Limit (mAbs)	3500
Cal. Reps Range (%)	***
Min Cal. Response (mAbs)	***
Cal. Factor Change (%)	***
M-Point Curve Fit (%)	***
Reagent Blank	ON
Auto R-Blank by Bottle	OFF

* Lot dependent
 ** operator definable
 *** optional
 N/A not applicable to this test

